the new



CONTOURA

by INMODE

YMNOON THOUSE

 $\overline{\mathsf{PLUS}}$

Full Body Contouring Experience



CONTOURA IS THE RESHAPING SOLUTION TO MAXIMIZE PATIENT RESULTS

Body contouring is one of the most popular aesthetic treatments. In many cases, patients are looking for improvements in body shape without surgical intervention. The Contoura is a non-invasive body contouring solution that will provide the ultimate in lunchtime procedures.

Contoura addresses:

»Body Reshaping

»Cellulite Reduction

»Skin Tightening

The Contoura is the ONLY non-invasive body contouring device proven in peer reviewed and published human studies to both permanently kill adipose tissue and contract the skin.

Provide your patients with clear and clinically proven results, so they can have the beautiful curves they want, from top to bottom.

BENEFITS TO YOUR PRACTICE

- » Contoura has zero consumable costs or pulse counters. Unlike other body contouring devices (which may carry high per procedure cost), Contoura is a business solution that allows you to grow your patient base without having to worry about ongoing fees.
- » The use of radio frequency allows owners to provide both skin contraction and fat treatments, of all body areas, in one device. This provides multiple revenue streams to a practice's bottom line.

Obesity trends are at an all-time high. According to the American

Heart Association, one in three US adults are obese. The average

American man has a 5 inch increase in waist size (now 40 inches) and the average American woman has a 7 inch increase in waist size (now 37 inches).

Contoura addresses this shift in the market.

GET THE ULTIMATE IN RESHAPING WITH A POWERFUL CONTOURING TOOL

| BODY f x | The first and only non-invasive body contouring technology proven in peer reviewed and published human studies to both permanently kill adipose tissue and contract the skin. |
|----------|---|
| MINIfx | Perfectly sized to address problematic fatty tissue in smaller areas like the arm, neck and knees. |
| PLUS | The first auto-adjusting thermal skin treatment for collagen remodeling and skin contraction all areas of the body. |

A.C.E. Technology

Acquire. Control. Extend. Technology to safely optimize peak energy

» ACQUIRE:

Smart RF thermistors built into the handpiece reads skin temperature 1000x per second to "ACQUIRE" therapeutic skin contraction temperatures without overheating. Measurements are clearly displayed on the system's control screen and audible indicators alert the operator to critical changes in the temperature and impedance profile in real time.

» CONTROL:

Unlike ANY other thermal skin enhancement technology to date, the A.C.E. "CONTROL" feature is a one-of-a-kind software algorithm that allows the clinician unprecedented safety and control by providing programmable "CUT-OFF" temperature between 35°-43°C.

» EXTEND:

Clinical evidence indicates that prolonged exposure at (41°C - 43°C) is advantageous for skin tightening. With the A.C.E. "EXTEND" feature you can now expose your patients skin, adipose tissue, or both to clinically therapeutic temperatures for as long as you want without the fear of adverse incidents in order to maximize your patients' outcomes and drastically reduce or eliminate non-responders.



BodyFX & MiniFX are a breakthrough third generation non-invasive body contouring technology for the treatment of the buttocks, abdomen, flanks and thighs. Not only does BodyFX selectively target problematic deep fatty tissue, but it also provides radio frequency heat to the surface of the skin, which aids the total area body contouring. With multiple peer reviewed publications, BodyFX is a proven long term solution with permanent fat cell death.

How does it work?

- Radio-Frequency provides precise and optimal heating of the skin for tightening and body shaping
- Vacuum massage applies gentle pressure allowing for maximal depth treatment
- A.C.E. Technology works in unison with the above to optimize ideal heating conditions for a safe and effective treatment



Apoptosis Study in Laser In Surgery Medicine

- Fat cell apoptosis: 30% of adipocytes, in a 2.5 cm treatment zone, stained positive for APAF-1, a validated marker for apoptosis
- Significant increase in collagen synthesis (neocollagenesis) representing an average increase of 13.7%



Plus is a non-invasive, sub-dermal heating applicator for collagen remodeling. RF current flows uniformly between the electrodes to provide a comfortable thermal experience. Built-in temperature controls allow for the attainment of clinically proven optimal temperatures, and the exclusive capacity for prolonged exposure at therapeutic temperatures.



$\mathbf{BODY} f x$

| Chamber Size: | 2 inches x 1 inch |
|----------------------|---|
| RF Output Frequency: | Up to 50 W |
| Output Frequency: | 1 MHz |
| Negative Pressure: | Up to 500 mBar |
| Temperature Cut-Off: | Real time temperature monitoring and cut-off. Operator selectable up to 43°C. |





DR. D. DUNCAN After 8 Treatments





DR. J. DIAMOND After 6 Treatments






DR. J. DIAMOND

After 6 Treatments





SERENITY CREEK MED SPA BodyFX & Forma



$\min f x$

| Chamber Size: | 1 inch x 1 inch |
|----------------------|---|
| RF Output Frequency: | Up to 50 W |
| Output Frequency: | 1 MHz |
| Negative Pressure: | Up to 500 mBar |
| Temperature Cut-Off: | Real time temperature monitoring and cut-off. Operator selectable up to 43°C. |









INMODE

DR. J. HELLMAN



PLUS

| RF Output Frequency: | Up to 65 W |
|----------------------|---|
| Output Frequency: | 1 MHz |
| Temperature Cut-Off: | Real time temperature measurements set by operator up to 43°C |





DR. B. RIDENOUR After 6 Treatments





DR. S. MULHOLLAND After 8 Treatments





DR. B. RIDENOUR After 6 Treatments





F. ORMONDE, MA After 8 Treatments







MINI fx



PLUS



OR CALL US AT {855} 411 2639 FOR A CONSULTATION

V0C0B030417